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Immunohistochemical Experiment Report of Frozen Section

I Experimental equipments and reagents

1.Experimental equipments

Name	Manufacturer	Model
Frozen slicer	Thermo	CRYOSTAR NX50
Slide	Wanwu	
Cover glass	Jiangsu Shitai Experimental Equipment Co., LTD	10212432C
Microwave oven	Galanz Microwave oven electrical appliances co., LTD	P70D20TL-P4
Decolorization shaker	Wanwu	TSY-B
Vortex mixer	Wanwu	MX-F
Palm centrifuge	Wanwu	D1008E
Pipette gun	Dragon	KE0003087/KA0056573
Tissue pencil	Wanwu	WG1066-1
Refrigerator	Qingdao Haier Co., LTD	BCD-192TGN
Microscope	Nikon	E100

2.Main experimental reagents

Reagent	Manufacturer	Model	Dilution Ratio
OCT embedding agent	Wanwu	G6059-110ML	
Anhydrous ethanol	China National Pharmaceutical Group Chemical Reagent Co., LTD	100092183	
Xylene	Shanghai Lingfeng Chemical Reagent Co., LTD	1330-20-7	
N-butanol	China National Pharmaceutical Group	100052190	
EDTA(PH9.0)antigen retrieval buffer	Chemical Reagent Co., LTD Wanwu	G1203	

Phosphate buffer saline	Wanwu	G0002
4% Paraformaldehyde	Wanwu	G1101
	China National	
3% Hydrogen peroxide	Pharmaceutical Group	10011208
	Chemical Reagent Co., LTD	
Bovine serum albumin	Wanwu	G5001
Normal rabbit serum	Wanwu	G1209
Hematoxylin stain solution	Wanwu	G1004
Hematoxylin differentiation solution	Wanwu	G1309
Hematoxylin returning blue solution	Wanwu	G1340
Neutral gum	Wanwu	G1403
Primary antibody:		
Secondary antibody:		
Immunohistochemical kit	Wanwu	G1211
DAB chromogenic agent		

II Immunohistochemical experimental procedure of frozen section

1. Frozen section fixation: the frozen sections are dried at room temperature, baked in an oven at 37°C for 10-20 minutes, fixed in 4% paraformaldehyde for 20 minutes, shaken on a decolorizing shaker in PBS(PH7.4) 3 times for 5 minutes each.
2. Antigen retrieval: the tissue sections are placed in a repair box filled with citric acid(PH9.0) antigen retrieval buffer for antigen retrieval in a microwave oven. Then, medium power for heating 5 minutes, power cut interval 5 minutes to low power for heating 5 minutes. During this process, excessive evaporation of buffer should be prevented and the sections should not be allowed to dry. After natural cooling, the sections are placed in PBS and shaken on the decolorization shaker 3 times for 5 minutes each. (Repair fluid and conditions are determined by tissue).
3. Blocking endogenous peroxidase activity: the sections are placed in 3% hydrogen peroxide and incubated at room temperature in darkness for 25 minutes. The sections are placed in PBS(PH7.4) and shaken on a decolorizing shaper 3 times for 5 minutes each.
4. Serum sealing: after the sections are slightly dried, a tissue pencil is used to draw a circle around the tissue (to prevent the antibodies from flowing away), and 3%BSA or 10% normal rabbit serum is added to cover the tissues evenly within the circle, and the tissues are sealed for 30 minutes at room temperature.(Primary antibody is sealed with 10% normal rabbit serum

- from goat source and other sources are blocked with 3%BSA).
5. Primary antibody incubation: the blocking solution is gently removed, the primary antibody prepared with PBS(PH7.4) in a certain proportion is added to the sections, and the sections are placed flat in a wet box and incubated overnight at 4°C. (Add a small amount of water in the wet box to prevent evaporation of antibodies).
 6. Secondary antibody incubation: the sections are placed in PBS(PH7.4) and washed by shaking on the decolorizing shaker 3 times for 5 minutes each. After the sections are slightly shaken and dried, the tissues are covered with secondary antibody (HRP labeled) from the corresponding species of primary antibody and incubated at room temperature for 50 minutes.
 7. DAB chromogenic reaction: the sections are placed in PBS(PH7.4) and shaken on the decoloring shaker 3 times for 5 minutes each. DAB color developing solution newly prepared is added in the circle after the sections are slightly dried. The color developing time is controlled under the microscope. The positive is brownish yellow. Rinse the sections with tap water to stop the reaction..
 8. Nucleus counterstaining: the sections are counterstained with hematoxylin stain solution for about 3 minutes; washed with tap water; differentiated with hematoxylin differentiation solution for several seconds; washed with tap water; treated with hematoxylin returning blue solution; washed with running water.
 9. Dehydration and mounting: place the section in 75% alcohol for 5 minutes--85% alcohol for 5 minutes--absolute ethanol I 5 minutes--anhydrous ethanol II 5 minutes--n-butanol 5 minutes--xylene I 5 minutes, dehydrated and transparent, remove the sections from xylene and let them dry slightly, then mount the sections with neutral gum.
 10. Visualize staining of tissue under a microscope, acquisitive and analysis image.

III Interpretation of the immunohistochemical results of frozen sections

The nucleus of hematoxylin stained is blue, and the positive expression of DAB is brownish yellow.