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## Chromatin Immunoprecipitation protocol

### 1、Experiment instruments

Name	Company	Model number
Microplate Reader	Rayto	RT-6100
Homogenizer	Wanwu	KZ-II
Centrifuge	DRAGONLAB	D3024R
Mini Centrifuge	Wanwu	D1008E
Vortex mixer	Wanwu	MX-F
Magnetic stirrers	Wanwu	MS-PB
Decolorizing shaker	Wanwu	TSY-B
Thermostatic Water Bath	Jiangyan Tianli Medical Devices Co., Ltd.	TL-420D
Camera obscura	Wanwu	WGA0017
Darkroom lamp	Wanwu	WGA0016
Photographic film	Wanwu	WGJP0001
Real-time qPCR	ABI	Stepone plus
Bechtop	AIRTECH	SW-CJ-1FD
Gel imaging system	Shanghai Tianneng Technology Co., LTD	Tanon-1600R
Ice Maker	sigma	SPR 80
Scanister	EPSON	V370
Gray analysis software	Alpha Innotech	alphaEaseFC
Image-analysis software	Adobe	Adobe PhotoShop
Ultrasonic Cell Disruptor	Ningbo Xinyi ultrasonic Equipment Co. LTD	JY 92-IIN

### 2、Reagents

Reagents	Company	Catlog
IP Lysis Buffer	Wanwu	G2038-100ML
Protein A/G-beads	Millipore	IP05
50*cocktail	Wanwu	G2006
PMSF (100mM)	Wanwu	G2008
Phosphoproteinase inhibitor	Wanwu	G2007
BCA protein quantitative detection kit	Wanwu	G2026
5*/2*loading buffer	Wanwu	G2013/G2031

---

SDS-PAGE kits	Wanwu	G2003
Protein Marker	Wanwu	26617
TRIS	Wanwu	G5009-50
Glycine	Wanwu	G5010-500G
SDS	Wanwu	G5037-100G
Milk	Wanwu	G5002-100G
TWEEN 20	Wanwu	WGT8220
ECL	Wanwu	G2014
Developer and fixer	Wanwu	G2019
PeroXidase-conjugated Goat Anti- Rabbit IgG(H+L)	Wanwu	GB23303
PeroXidase-conjugated Donkey Anti-Goat IgG(H+L)	Wanwu	GB23404
PeroXidase-conjugated Goat Anti-Mouse IgG(H+L)	Wanwu	GB23301
PeroXidase-conjugated Goat Anti-Rat IgG(H+L)	Wanwu	GB23302
Transfer buffer	Wanwu	G2017
Electrophoresis buffer	Wanwu	G2018
TBS buffer	Wanwu	G0001
Agarose	Sinopharm Group Chemical	G5056
Anhydrous ethanol	Reagent Co. LTD	10009218
2×SYBR Green qPCR Master Mix (High ROX)	Wanwu	G3322

---

### 3、 Protocol For ChIP

#### Day 1:

##### 3.1 Formaldehyde cross-linking and ultrasonic disruption of cells

3.1.1. For cells: measure the volume of the medium, and add 37% formaldehyde until a final concentration of 1% is reached.

3.1.2. Incubate for 10 minutes at 37 °C

3.1.3. Terminate the cross-linking: Block the reaction with 500 µl Glycine 2.5 M (final concentration 125 mM). Incubate for 5 minutes at room temperature.

3.1.4. Discard the medium and wash the cells 3 times with ice cold PBS;

3.1.5. Using a cell scraper scrape off the cells in PBS, centrifuge at 2000g/min for 5min, and then

discard the supernatant;

3.1.6. Breaking the cells in IP Lysis Buffer supplemented with protease inhibitor(the volume of IP Lysate Buffer depends on the amount of cells).

3.1.7. To facilitate the cell membrane breaking, the cells were repeatedly blown with a pipette (or shaken on a vortex mixer);

3.1.8. Ultrasonic Crushing: Ningbo Xinzhi JY 92-IIN, 2# probe, 30% power, 4.5s ultrasound, 10.5s interval, 5min, and 20 times. During the ultrasonic process, please be careful to keep the sample in the ice bath, avoid ultrasonic bubbles .

### **3.2 Re-clear samples, and IP overnight**

3.2.1. After ultrasonic breaking, centrifuge at 12000rpm, 4°C for 10min. Remove supernatant into a fresh tube. Transfer 90μL lysis solution for Input, store the others at -80 °C.

3.2.2. Take 40 μ l ultrasonic crushing product as input, add 10 μ L 5 \* loading buffer, heat-denature it and perform WB detection to confirm that the sample contains the target protein. The remaining 50 μ l product was added with 2.5 μ l protease K and 2 μL 5M NaCl (the final concentration of NaCl was 0.2M), and the solution was decrosslinked at 55°C overnight;

3.2.3. Measure the nucleic acid concentration .Take part of the sample for PCR amplification, and then run agarose electrophoresis to detect the effect of ultrasonic , and confirm whether the sample contains the target DNA;

3.2.4. After confirming the Input result, take 100 μl of the ultrasonically crushed product , add 900 μl of ChIP Dilution Buffer containing 1 mM-PMSF and 20 μl of 50×PIC (cocktail). Add 60 μl of Protein A+G Agarose/Salmon Sperm DNA for each sample. Incubate for 1 hours at 4°C in a rotating wheel.

3.2.5. Incubate for 10 minutes at 4°C , centrifuge at 4000 rpm for 5 minutes;

3.2.6. Transfer the sample into two 1.5mL EP tubes.For the IgG tubes,add 1μg of IgG (of the same species) . For the IPs, use 2 μg of antibody per tube. Incubate overnight at 4°C in a rotating wheel.

### **Day 2:**

#### **3.3 Wash samples, elute, and reverse crosslink.**

3.3.1. Add 200μl Protein A+G Agarose/Salmon Sperm DNA to each tubes . Rotate for 2h at 4 °C to pull down immuno complexes.

3.3.2. Let all tubes standing at 4°C for 10min. Centrifuge at 4000 rpm for 1 minutes. Discard the supernatant.

3.3.3. Wash the precipitate complex with the following solution in turn. Cleaning steps: add the solution, turn it upside down slightly, centrifuge at 4000rpm for 3min, and remove the supernatant.

Washing solution

1x low salt immune complex wash buffer.

1x high salt immune complex wash buffer.

1x LiCl immune complex wash buffer.

2x TE buffer.

The samples can be used for Western Blot analysis. Add 30 $\mu$ l 1X loading buffer to 30 $\mu$ l beads, then denature them for WB.

3.3.4. Freshly prepare elution buffer (0.5g SDS + 0.42g NaHCO<sub>3</sub> + 50ml Purified Water). Add 200  $\mu$ L of elution buffer to elute protein;

3.3.5. Add 20 $\mu$ l 5 M NaCl and 10 $\mu$ l proteinase K to the eluents. And heat all samples at 55 °C overnight to reverse the crosslinking.

### **Day 3:**

#### **3.4 Purify DNA with Spin Column**

3.4.1. Add 200  $\mu$ l Buffer GB to the sample, mix thoroughly by vortex, and incubate at 70°C for 10min to yield a homogeneous solution. Briefly centrifuge the 1.5 ml microcentrifuge tube to remove drops from the inside of the tube.

3.4.2. Add 200 $\mu$ l ethanol (96-100%) to the sample, and mix thoroughly by vortex for 15 s. Briefly centrifuge the 1.5 ml microcentrifuge tube to remove drops from the inside of the tube.

3.4.3. Transfer the mixture from step 2 into the Spin Column CB3 (in a 2 ml collection tube) and centrifuge at 12,000 rpm (~13,400  $\times$  g) for 30s. Discard flow-through and place the spin column into the collection tube.

3.4.4. Add 500  $\mu$ l Buffer GD (Ensure ethanol (96-100%) has been added) to Spin Column CB3, and centrifuge at 12,000 rpm (~13,400  $\times$  g) for 30s, then discard the flow-through and place the spin column into the collection tube

3.4.5. Add 700 $\mu$ l Buffer PW (Ensure ethanol (96-100%) has been added) to Spin Column CB3, and centrifuge at 12,000 rpm (~13,400  $\times$  g) for 30s. Discard the flow-through and place the spin column into the collection tube.

3.4.6. Add 500 $\mu$ l Buffer PW to Spin Column CB3, and centrifuge at 12,000 rpm (~13,400  $\times$  g) for 30s. Discard the flow-through and place the spin column into the collection tube.

3.4.7. Centrifuge at 12,000 rpm (~13,400  $\times$  g) for 2 min, and place it at room temperature for 10min to dry the membrane completely.

3.4.8. Transfer the Spin Column CB3 in a new clean 1.5 ml microcentrifuge tube, and add 50-200  $\mu$ l Buffer TE directly to the center of the membrane. Incubate at room temperature (15-25°C) for 2min, and then centrifuge for 2min at 12,000 rpm (~13,400  $\times$  g)

#### **3.5 DNA analysis (Real-Time PCR or DNA Sequencing)**

##### **3.5.1. Real-Time PCR**

Using JASPAR or other database to predict the transcription factor binding sites. Design the primers according to the sites.

. Make a mixture in 200 $\mu$ l microcentrifuge tube

2 $\times$  qPCR Mix      10 $\mu$ l

---

2.5 $\mu$ M primer	2 $\mu$ l
Template	2 $\mu$ l
ddH <sub>2</sub> O	6 $\mu$ l

Run reaction

Pre-incubation	95°C, 10min
Amplification: 40 cycles of	95°C, 15s→60°C, 60s
Melting curve	60°C→95°C, the temperature grows 0.3°C every 15s

.Run PCR reaction on a 2.5% TAE gel and document.

### 3.5.2.DNA Sequencing

Collect DNA, measuring the concentration, and then send to a sequencing company.